


# Guide to Cell Modification & Re-Engineering Solutions



**HARVARD**  
APPARATUS

---

[www.harvardapparatus.com](http://www.harvardapparatus.com)

# Table of Contents

---

<b>Cell Diagram</b> .....	<b>2</b>
<b>Introduction</b> .....	<b>3</b>
<b>Select Proper Technique Table</b> .....	<b>4</b>
<b>Technique Definitions</b> .....	<b>5-6</b>
<b>Pneumatic Injectors</b>	
PLI-100 .....	7
PLI-90 .....	7
<b>Iontophoresis</b>	
Neuro Phore BH-2 System .....	8
<b>Mechanical Injector</b>	
Pump 11 Pico Plus .....	9
<b>Electroporation</b> .....	<b>10</b>
<b>Liposomes</b>	
Liposomat .....	11
<b>Reference</b> .....	<b>12</b>
<b>Harvard Apparatus Information</b> .....	<b>13</b>
<b>Technical Support</b> .....	<b>14</b>

# Cell Diagram

---

No matter what your cell biology application  
Harvard Apparatus has the solution


ENDOCYTOSIS

EXOCYTOSIS

INJECTION

TRANSFECTION

IMAGING



# Introduction

---

## The Harvard Apparatus Advantage...Serving cell biologists with exceptional products and support.

Introducing materials into cells such as, bacteria, DNA, drugs and ions, for the purpose of:

- Monitoring the cells behavioral to the materials
- Modifying and re-engineering the cells structure and function requires many different techniques, to get the cells to incorporate the materials into the cell to get desired effects

When you choose a techniques you must ask the following questions to select the optimal technique for your applications:

1. How many cells do I need for my experiment?
2. In what period of time do I need this number of cells?
3. Are my cells sensitive to certain materials in a negative way?
4. Of the cells modified do I need a high or low yield rate?
5. Are my cells temperature sensitive?
6. Are you going to keep them observed for longer than 4 hours?
7. Are my cells free floating or fixed?
8. How big is the substance going into the cell?
9. Do you only want some cells in a mixture handled?
10. Do I work with organs, tissue or cells clusters or individual cells?


# Select Proper Technique

Harvard Apparatus has many ways for you to handle the introduction of samples into the Cells that best meets your experimental goals. The following chart will assist you in your selection:


	Electroporation	Iontophoresis	Liposomes	Pneumatic Injectors	Mechanical
<b>Example</b>	<b>BTX</b>	<b>BH-2</b>	<b>Liposomat</b>	<b>PLI-100</b>	<b>Pico plus</b>
<b>Technical driving principal</b>	High Voltage Opens Pores	Charge thru Existing Pores	Fat Ball thru Cell Wall	Push in w/Air	Push in w/Liquid
Average Price	\$ 8000	\$ 10000	\$ 3000	\$ 4000	\$ 2000
One Cell at a Time	Not Typical	Not Typical	Not Typical	Yes	Yes
Transfect Hundreds of Cells Simultaneously	Yes	Yes	Yes	No	No
Transfect Thousands of Cells Simultaneously	Yes	No	Yes	No	No
Cell Container	Cuvette or In Vivo	In Vivo	Culture or In Vivo	Culture	Culture
Typical yields	30-70%	30-70%	50-90%	>90%	>90%
Gentle to Cells	Harsh to Low	Medium to Low	Low	Low	Low
Labor Intensive	Yes	Yes	Yes	No	No
Require Lots of Technique	Little	Little	Some	Yes	Yes
Need to Hold Cells Rigidly	No	No	No	Recommended	Recommended
Needs Cell Rigidly Held	Not Need	Not Need	Not Need	Yes	Yes
<b>Equipment:</b>					
Microscope	No	No	No	Yes	Yes
Capillary Glass Puller	No	No	No	Yes	Yes
Microforge	No	No	No	Yes	Yes
<b>Applications:</b>					
Large Particles and Molecules	Yes	No	No	Yes	Yes
Small Particles and Molecules	Yes	Yes	Yes	Yes	Yes
Bacteria	Yes	No	No	Yes	Yes
DNA	Yes	No	No	Yes	Yes
Virus	Yes	No	Yes	Yes	Yes
RNA	Yes	No	No	Yes	Yes
Drugs	Yes	Yes	Yes	Yes	Yes
Proteins	Yes	No	No	Yes	Yes
Plasmids	Yes	Yes	Yes	Yes	Yes

# Technique Definitions

**Iontophoresis** is a non-invasive method of propelling high concentrations of a charged substance, normally medication or bioactive-agents, transdermally by repulsive electromotive force using a small electrical charge applied to an iontophoretic chamber containing a similarly charged active agent and its vehicle. To clarify, one or two chambers are filled with a solution containing an active ingredient and its solvent, termed the vehicle. The positively charged chamber, termed the anode will repel a positively charged chemical, whilst the negatively charged chamber, termed the cathode, will repel a negatively charged chemical into the skin.




**Electroporation**, or electroporation, is a significant increase in the electrical conductivity and permeability of the cell plasma membrane caused by an externally applied electrical field. It is usually used in molecular biology as a way of introducing some substance into a cell, such as loading it with a molecular probe, a drug that can change the cell's function, or a piece of coding DNA. Pores are formed when the voltage across a plasma membrane exceeds its dielectric strength. If the strength of the applied electrical field and/or duration of exposure to it are properly chosen, the pores formed by the electrical pulse reseal after a short period of time, during which extracellular compounds have a chance to enter into the cell. However, excessive exposure of live cells to electrical fields can cause apoptosis and/or necrosis -the processes that result in cell death.




# Technique Definitions (cont.)

A **liposome** is a spherical vesicle with a membrane composed of a phospholipid and cholesterol bilayer (see on the right). Liposomes can be composed of naturally-derived phospholipids with mixed lipid chains (like egg phosphatidylethanolamine), or of pure surfactant components like DOPE (dioleoylphosphatidylethanolamine). Liposomes are used for drug delivery due to their unique properties. A



liposome encapsulates a region on aqueous solution inside a hydrophobic membrane; dissolved hydrophilic solutes can not readily pass through the lipids. Hydrophobic chemicals can be dissolved into the membrane, and in this way liposome can carry both hydrophobic molecules and hydrophilic molecules. To deliver the molecules to sites of action, the lipid bilayer can fuse with other bilayers such as the cell membrane, thus delivering the liposome contents. By making liposomes in a solution of DNA or drugs (which would normally be unable to diffuse through the membrane) they can be (indiscriminately) delivered past.

An **injector** is a pump-like device, which charges or discharges containers under pressure with suitable arrangements. Motive force is gained at the inlet from a



suitable gas or liquid that is under pressure. Pneumatic injectors can deliver 10–15L femtolitre of materials and specialized mechanical pumps are capable of 10–12L picolitre delivery volumes.

# Pneumatic Injectors

---

## PLI-100 and PLI-90 Pico-Injector Series

- 5 pressures: inject, balance, clear, fill and hold
- Popular applications:
  - Injection of mouse, frog, zebrafish and other oocytes
  - Extracellular brain injections
  - Injection of DNA, mRNA, microbeads, neurotransmitters, kinases and other proteins
- Most Published Injector

The PLI-100 Pico-Injector reliably delivers a wide range of volumes through micropipettes by applying a regulated pressure for a digitally set period of time. Compressed gas allows the user to deliver desired volumes from femtoliters to microliters while simultaneously holding a cell. Whether you need to do large injections into capillaries or very small injections into mammalian nuclei, the PLI-100 is well suited for your experiment.

The **PLI-100** has become a favorite of prestigious national microinjection workshops like Cold Spring Harbor Laboratories and other researchers worldwide. Other companies have tried to design similar systems, but the PLI-100 remains unparalleled in terms of ease of use, durability, precision, and cost. Pico-Injector, like the reliable PLI-100 Pico-Injector, controls the precise and reproducible regulation of injection pressure and time.


Because the **PLI-90** features only the injection, balance, and clearing positive pressures, it is a lower cost alternative to the PLI-100. This simplicity makes it even easier to use than the PLI-100.



# Iontophoresis

## Neuro Phore BH-2

- **Modular design; buy what you need now and add additional modules later as your protocol evolves**
- **Accommodates 7 barrel micropipettes supports the most complex and demanding injection and recording protocols**
- **Extremely low noise,  $\pm 105$  volts compliance, allows rapid iontophoretic injection with even the smallest micropipette tips**
- **Automatic current neutralization for minimal electrical artifacts**
- **Digitally controlled eject and pause timing utilizes easy to read and set digital panel switches on each iontophoretic or pressure module/channel**



The **Neuro Phore BH-2 System** was developed under the guidance of active researchers with extensive experience in iontophoresis techniques. These researchers needed a system to provide precise stimulation and quantitative control for ejection of drugs in their pharmacological studies of drug evoked responses such as neuro-synaptic discharges, contraction, and changes in chemical concentration. What emerged was a reliable, accurate, easy to use, Neuro Phore BH-2 system that is capable of accommodating high impedance multi-barreled micro-electrode pipettes.

# Mechanical Injector

---

- Infuse/ withdraw capability
- Dual syringes for broad flow rate range
- Easy-to-read display, the injection rate can be changed from 0.0550  $\mu\text{l}/\text{sec}$  to 0.0073  $\text{ml}/\text{sec}$
- Remote to micromanipulator injections are easy non-obstructed viewing and no heavy weight to hinder positioning
- This pump has exceptional milliliter, microliter, nanoliter and picoliter smooth and accurate flow, easily set through the keyboard
- The injection can be controlled by initiating the manual start button, RS232 through your computer or with the included foot pedal
- Operating parameters are set with the membrane keypad and VFD display
- From the keypad the user can set pump to:
  - Infusion mode
  - Withdrawal mode (reversing switch on back panel)
- All micro tubing and connectors are available
- Volume mode -enter a target volume, pump will stop when value achieved
- Adjustable infusion limit switch and adjustable withdraw mechanical stop



The **Pump 11 Pico Plus** is the lowest flow rate pump that Harvard Apparatus manufactures. It can accurately delivery volumes down to 3.3 picoliters/minute. It is ideal for applications such as cellular injections or for providing the micro flow required for FIA or capillary LC. It has both infusion and withdrawal capabilities as well as RS-232 communications.

## **Applications:**

- Cellular injection i.e. oocytes
- Micro flow for FIA or capillary LC
- Micro reaction delivery
- HPLC post-column reaction pump

# Electroporation

- Specialty electrodes
- Cuvettes
- High throughput  
96-well plates

Electroporation is a techniques used to perform research on embryonic stem cells (HESC), adult stem cells, somatic cell nuclear transfer, therapeutic cloning and cancer research.

## Applications:

- Transfection of Mammalian Cells
- In Vitro, In Vivo, Ex Vivo and In Ovo

## Experiments

- Nuclear Transfer
- Plant Tissue and Protoplast Transfection
- Transformation of Bacteria and Yeast



# Liposomes


---

- **Liposomes are used for studying protein-liposome interactions with compounds such as serum lipoproteins, lectins, toxins and clotting components.**
- **Entrapment of anti-tumor agents, anti-microbial drugs, anti-inflammatory and immunomodulatory agents and CNS-active drugs. The therapeutic applications are presently include cancer therapy, arthritis, metal chelation therapy, enzyme replacement therapy, hemophilia (factor VIII), myocardial infarction, and as radiopharmaceutical markers.**
- **Liposomes are used for studying in vitro and in vivo liposome-cell interactions.**
- **They are used for reconstitution experiments especially for ion transport systems.**
- **In molecular biology, liposomes are used for the mediated delivery of macromolecules into eukaryotic cells.**
- **The use of liposomes in gene-based technology and in broader biotechnology and pharmaceutical applications is increasing day by day.**

Today liposomes are an important part of biological, pharmaceutical, and medical research. Because liposomes are the most effective carriers for the introduction of many different types of agents into cells, the applications of liposome-based samples and products are extremely wide.

## Liposomat

- **The Liposomat is ideal for the preparation of liposomes of volumes from 3 ml to 50 ml or higher.**
- **The system has two serpentine channels superimposed on each other and separated by a membrane.**
- **Each channel has a volume of 3 ml and a length of 3 meters. The mixed lipid/detergent micelles run through one of the channels while the buffer flows through the other channel.**
- **Due to controlled dialysis and the high surface area in the system, liposomes can be formed within 30 minutes.**



# Reference

---

## Reference Articles:

Measurement of volume injected into individual cells by quantitative fluorescence microscopy,  
Lee GM.,  
J Cell Sci. 1989 Nov;94 ( Pt 3):443-7

Information Fusing Control of the DNA Microinjection Volume Based on Stochastic Fuzzy Neural  
Network,  
Zhang Ling, Yu Yongquan, ZengBi  
2006-02-21 09:22:58.0

Fluid Dynamics in Conically Tapered Microneedles,  
WijayaMartanto, Stephen M. Baisch, Elizabeth A. Costner, Mark R. Prausnitz, Marc K. Smith  
AIChEJournal, Vol.51, No.6, 1599-1607, 2005

Controlled micro release of pharmacological agents: measurements of volume ejected in vitro through  
fine tipped glass microelectrodes by pressure,  
M. SAKAI\*, BARBARA E. SWARTZ and C. D. WOODY  
Neuropharmacology. 1979 Feb;18(2):209-13.

Easy-to-use equipment for the accurate microinjection of nanoliter volumes into the nuclei of  
amphibian oocytes,  
Stephens DL, Miller TJ, Silver L, ZipserD, Mertz JE.  
Anal Biochem. 1981 Jul 1;114(2):299-309.

# Harvard Apparatus Information

---

## GETTING TO KNOW HARVARD APPARATUS!

Since 1901, **Harvard Apparatus** has been serving the global bioresearch community by developing solutions for the thousands of experimental challenges in animal, organ, tissue, cellular, and molecular bioresearch. Harvard Apparatus offers "Total Solutions" which are utilized in the area of studies of: Physiology, Electrophysiology, Biology, Neurology, Cardiology, Toxicology, Pharmacology, Molecular Biology. Harvard Apparatus is made up of seven groups developing specialized products:

### WHOLE ANIMAL

Harvard Apparatus

Syringe pumps for small animal infusions and MS calibration; Small animal ventilators and anesthesia set-ups

Harvard Apparatus LTD

Harvard Apparatus/Panlab

Kymographs, biosensors and plethysmographs  
Behavioral research Systems

### MOLECULAR BIOLOGY

Harvard Apparatus

Sample Preparation -

CPK molecule models to SPE columns for de-salting and equilibrium dialysis in 1 or 96 well plates for T3/T4 measurements or protein-lig and binding studies

### ORGANS AND TISSUES

Harvard Apparatus/

Hugo Sachs -

Tools for toxicological and pharmacokinetic studies utilizing isolated organ systems and perfusion baths, data acquisition systems, amplifiers and transducers

### CELLS

Harvard Apparatus/

Warner Instruments

Observing, modifying, recording and maintaining cells under a microscope or ion transport studies through a cell membrane i.e. Ussing Chambers

Harvard Apparatus/BTX

A full family of electroporation and electrofusion Products Our applications experts look forward to assisting you at [www.harvardapparatus.com](http://www.harvardapparatus.com), click on contact us to find local support.

Ask for any or all of our free catalogs: Cell Biology; Physiology, Pharmacology, Toxicology; Sample Preparation, Infusion Pumps, Behavioral Sciences, Ventilation /Anesthesia, and Surgical Supplies and Equipments.

# Technical Support

---

## Technical support for pre-sales assistance and post sales support

### Headquarters

#### UNITED STATES



#### Harvard Apparatus

Attn: Customer Service  
84 October Hill Road  
Holliston, Massachusetts 01746, USA

phone **508.893.8999**

fax **508.429.5732**

e-mail **bioscience@harvardapparatus.com**

website **www.harvardapparatus.com**

Harvard Apparatus believe's in providing you the best technical support. Our support specialist can assist you with pre-sales recommendations and assist you In configuring complex infusion or anesthesia set-ups; selecting the appropriate surgical equipment, or selecting the best pump for microdialysis. We assist in making sure your purchase performs your tasks to your satisfaction

### Subsidiaries

#### Canada



#### Harvard Apparatus Canada

Attn: Sales Department  
6010 Vanden Abeele  
Saint-Laurent, Quebec H4S 1R9, Canada

phone **514.335.0792**

**800.361.1905** (canada only)

fax **514.335.3482**

e-mail **sales@harvardapparatus.ca**

website **www.harvardapparatus.ca**

#### SPAIN



#### Panlab, S.L.

#### Harvard Apparatus Spain

C/Energia, 112  
08940 Cornellà, Barcelona  
Spain

phone **34.934.750.697** (International Sales)

phone **934.190.709** (Sales in Spain)

fax **34.934.750.699**

e-mail **info@panlab.com**

website **www.panlab.com**

#### Germany



#### Hugo Sachs Elektronik

#### Harvard Apparatus, GmbH

Gruenstrasse 1  
D-79232 March-Hugstetten, Germany

phone **49.7665.92000**

fax **49.7665.920090**

e-mail **info@hugo-sachs.de**

website **www.hugo-sachs.de**

#### United Kingdom



#### Harvard Apparatus, Ltd.

Attn: Sales Department, Fircroft Way,  
Edenbridge, Kent TN8 6HE, United Kingdom

phone **44.1732.864001**

fax **44.1732.863356**

e-mail **sales@harvardapparatus.co.uk**

website **www.harvardapparatus.co.uk**

#### France



#### Harvard Apparatus, S.A.R.L.

Attn: Sales Department  
6 Avenue des Andes, Miniparc – Bat. 8  
91952 Les Ulis Cedex, France

phone **33.1.64.46.00.85**

Fax **33.1.64.46.94.38**

e-mail **info@harvardapparatus.fr**

website **www.harvardapparatus.fr**