




Document Title: Order form for custom made probes & guides

Customer:	Order No (to be filled in by CMA):
No of pkg	Customer's order No/Ref:

Choose the type of probe and fill in the missing information in the schematic drawing below. Minimum order qty of probes is 3 pkg. For guides see page 3. For detailed information regarding materials, etc contact CMA Microdialysis AB.

CMA 7 (3 pcs/pkg) <input type="checkbox"/> Ref. No. 8010391 3 pkg <input type="checkbox"/> Ref. No. 8010392 4-10 pkg <input type="checkbox"/> Ref. No. 8010393 >10 pkg	CMA 8 (3 pcs/pkg) <input type="checkbox"/> Ref. No. 8012401 3 pkg <input type="checkbox"/> Ref. No. 8012402 4-10 pkg <input type="checkbox"/> Ref. No. 8012403 >10 pkg
<p> <input type="checkbox"/> normal 200 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> normal 200 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> normal 7 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> Metal free </p> <p> <i>*For longer membrane than 2mm, use CMA 11</i> <i>** Metal free option not available</i> </p> <p> Membrane <input type="checkbox"/> 6k Da CO CUPR <input type="checkbox"/> 55k Da CO PES* <input type="checkbox"/> 500k Da CO PES** <input type="checkbox"/> 2M Da CO PES** </p>	<p> <input type="checkbox"/> normal 200 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> normal 200 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> normal 7 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> Metal free </p> <p> membrane <input type="checkbox"/> 20k Da CO PAES <input type="checkbox"/> 100k Da CO PES </p>
CMA 11 (3 pcs/pkg) <input type="checkbox"/> Ref. No. 8010394 3 pkg <input type="checkbox"/> Ref. No. 8010395 4-10 pkg <input type="checkbox"/> Ref. No. 8010396 >10 pkg	CMA 12 (3 pcs/pkg) <input type="checkbox"/> Ref. No. 8010397 3 pkg <input type="checkbox"/> Ref. No. 8010398 4-10 pkg <input type="checkbox"/> Ref. No. 8010399 >10 pkg
<p> <input type="checkbox"/> normal 14 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> Metal free </p> <p> <i>** Metal free option not available</i> </p> <p> membrane <input type="checkbox"/> 6k Da CO CUPR <input type="checkbox"/> 55k Da CO PES <input type="checkbox"/> 500k Da CO PES** <input type="checkbox"/> 2M Da CO PES** </p>	<p> <input type="checkbox"/> normal 14 mm <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> _____ mm </p> <p> <input type="checkbox"/> Metal free </p> <p> membrane <input type="checkbox"/> 20k Da CO PAES <input type="checkbox"/> 100k Da CO PES </p>

Document Title: Order form for custom made probes & guides

<p style="text-align: center;">CMA 20 (3 pcs/pkg)</p> <p><input type="checkbox"/> Ref. No. 8010400 3 pkg <input type="checkbox"/> Ref. No. 8010401 4-10 pkg <input type="checkbox"/> Ref. No. 8010402 >10 pkg</p> <div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> <p>outlet length <input type="checkbox"/> normal 200 mm <input type="checkbox"/> _____ mm</p> </div> <div style="width: 45%;"> <p>inlet length <input type="checkbox"/> normal 200 mm <input type="checkbox"/> _____ mm</p> </div> </div>  <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"> <p>shaft length <input type="checkbox"/> _____ mm</p> </div> <div style="width: 45%;"> <p>membrane length <input type="checkbox"/> _____ mm</p> </div> </div> <div style="margin-top: 10px;"> <p>membrane <input type="checkbox"/> 20k Da CO PAES <input type="checkbox"/> 100k Da CO PES</p> </div> <div style="border: 1px solid black; padding: 5px; margin-top: 10px; width: fit-content; margin-left: auto;"> <p>Normal length of shaft/membrane <input type="checkbox"/> 14 / 10 mm <input type="checkbox"/> 20 / 4 mm</p> </div>	<p style="text-align: center;">CMA 30 (4 pcs/pkg)</p> <p><input type="checkbox"/> Ref. No. 8010498 3 pkg <input type="checkbox"/> Ref. No. 8010499 4-10 pkg <input type="checkbox"/> Ref. No. 8010500 >10 pkg</p> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"> <p>Syringe connector <input type="checkbox"/> connector <input type="checkbox"/> no connector</p> </div> <div style="width: 45%;"> <p>membrane length <input type="checkbox"/> _____ mm</p> </div> </div>  <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"> <p>inlet length <input type="checkbox"/> normal 245 mm <input type="checkbox"/> _____ mm</p> </div> <div style="width: 45%;"> <p>outlet length <input type="checkbox"/> normal 245 mm <input type="checkbox"/> _____ mm</p> </div> </div> <div style="margin-top: 10px;"> <p>membrane <input type="checkbox"/> 6k Da CO CUPR</p> </div> <div style="border: 1px solid black; padding: 5px; margin-top: 10px; width: fit-content; margin-left: auto;"> <p>Maximal total length with inlet/outlet: 500 mm. Length of membrane: 4 or 10 mm.</p> </div>
<p style="text-align: center;">CMA 31(4 pcs/pkg)</p> <p><input type="checkbox"/> Ref. No. 8012406 3 pkg <input type="checkbox"/> Ref. No. 8012407 4-10 pkg <input type="checkbox"/> Ref. No. 8012408 >10 pkg</p> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"> <p>Syringe connector <input type="checkbox"/> connector <input type="checkbox"/> no connector</p> </div> <div style="width: 45%;"> <p>membrane length <input type="checkbox"/> 10 mm</p> </div> </div>  <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"> <p>inlet length <input type="checkbox"/> normal 350 mm <input type="checkbox"/> _____ mm</p> </div> <div style="width: 45%;"> <p>outlet length <input type="checkbox"/> normal 100 mm <input type="checkbox"/> _____ mm</p> </div> </div> <div style="margin-top: 10px;"> <p>membrane <input type="checkbox"/> 55k Da CO PES <input type="checkbox"/> 500k Da CO PES <input type="checkbox"/> 2M Da CO PES</p> </div> <div style="border: 1px solid black; padding: 5px; margin-top: 10px; width: fit-content; margin-left: auto;"> <p>Maximal total length with inlet/outlet: 500 mm. Length of membrane: 10 mm.</p> </div>	<p>Additional notes:</p>

Document Title: Order form for custom made probes & guides

Guides

CMA 7 guide

Guide shaft will be 2 mm shorter than probe shaft. The minimum shaft length is 3 mm, suitable to a probe with shaft 5 mm.*


CMA 8 guide

Guide shaft will be 2 mm shorter than probe shaft. The minimum shaft length is 5 mm, suitable to a probe with shaft 7 mm.*

CMA 11 & 12 guide

Guide shaft will be 4 mm shorter than probe shaft. The minimum shaft length is 3 mm, suitable to a probe with shaft 7 mm.*

* The length of the Guide Shaft determines the maximum implantation depth for the start of the probe membrane



CMA 7 guide

Ref. No. 8010313 3/pkg

Ref. No. 8010314 30/pkg

Specify details below in "additional notes"

CMA 8 guide

Ref. No. 8012404 3/pkg

Ref. No. 8012405 30/pkg

Specify details below in "additional notes"

CMA 11 guide

Ref. No. 8309029 3/pkg

Ref. No. 8309030 30/pkg

Specify details below in "additional notes"

CMA 12 guide

Ref. No.8309008 3/pkg

Ref. No. 8309009 30/pkg

Specify details below in "additional notes"

Additional notes:

For CMA Microdialysis production: The following special instructions shall be used in the production of the above probes:

Production manager, sign/date: